

Product Datasheet

FCGR2A Antibody (APC) (orb1272053)

Description	FCGR2A Antibody (APC)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Target	FCGR2A
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P12318
Expiration Date	12 months from date of receipt.