


Product Datasheet

ITGAX Antibody (PE/Cy7) (orb1272097)

Description	ITGAX Antibody (PE/Cy7)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Target	ITGAX
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P20702
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...