

Product Datasheet

Itgax Antibody (PE/Cy7) (orb1272137)

Description	Itgax Antibody (PE/Cy7)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Target	Itgax
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	Q9QXH4
Expiration Date	12 months from date of receipt.