


Product Datasheet

ITGAM Antibody (FITC) (orb1272157)

| | |
|----------------------------|--|
| Description | ITGAM Antibody (FITC) |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | FITC |
| Tested Applications | FC |
| Target | ITGAM |
| Preservatives | Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2. |
| Form/Appearance | liquid |
| Concentration | batch dependent |
| Storage | The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | Mouse IgG1, kappa |
| Clonality | Monoclonal |
| Uniprot ID | P11215 |
| Expiration Date | 12 months from date of receipt. |



Human peripheral blood monocytes were st...