

Product Datasheet

Emr1 Antibody (Biotin) (orb1272189)

| | |
|----------------------------|--|
| Description | Emr1 Antibody (Biotin) |
| Species/Host | Rat |
| Reactivity | Mouse |
| Conjugation | Biotin |
| Tested Applications | FC |
| Target | Emr1 |
| Preservatives | Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2. |
| Form/Appearance | liquid |
| Concentration | batch dependent |
| Storage | The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | Rat IgG2a, kappa |
| Clonality | Monoclonal |
| Uniprot ID | Q61549 |
| Expiration Date | 12 months from date of receipt. |