
Product Datasheet

CCL20 Antibody (Biotin) (orb1272362)

Description

CCL20 Antibody (Biotin)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Biotin

Tested

ELISA, WB

Applications
Immunogen

Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hMIP-3-alpha (human MIP-3-alpha).

Target

CCL20

Form/Appearance

Lyophilized

Concentration

batch dependent

Storage

MIP-3a antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes


ELISA: Sandwich: To detect hMIP-3α by sandwich ELISA (using 100 μL/well antibody solution) a concentration of 0.25 - 1.0 μg/mL of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our Polyclonal Anti-Human MIP-3α (XP-5240) as a capture antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hMIP-3α. **Western Blot:** To detect hMIP-3-α by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 μg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hMIP-3-α is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.

Clonality

Polyclonal

Uniprot ID
[P78556](#)
NCBI
[P78556](#)
Dilution Range

ELISA: Sandwich: To detect hMIP-3α by sandwich ELISA (using 100 μL/well antibody solution) a concentration of 0.25 - 1.0 μg/mL of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our Polyclonal Anti-Human MIP-3α (XP-5240) as a capture antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hMIP-3α. **Western Blot:** To detect



To detect Human MIP-3-alpha by sandwich ...