

Product Datasheet

NAMPT Antibody (orb1272636)

Description

NAMPT Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Produced from sera of rabbits immunized with highly pure recombinant Human Visfatin. Human Visfatin specific antibody was purified by affinity chromatography employing an immobilized Human Visfatin matrix.

Target

NAMPT

Form/Appearance

Lyophilized

Concentration

batch dependent

Storage

The lyophilized antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

ELISA: Sandwich: To detect Human Visfatin by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with Biorbyt's Biotinylated Anti-Human Visfatin (38-282) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human Visfatin. Western Blot: To detect Human Visfatin by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant Human Visfatin is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.

Clonality

Polyclonal

Uniprot ID


[P43490](#)

NCBI

[P43490](#)

Expiration Date


12 months from date of receipt.



To detect Human Visfatin by sandwich ELISA



To detect Human Visfatin by Western Blot



To detect Human Visfatin by Western Blot