
Product Datasheet

TNFSF18 Antibody (orb1272727)

Description	TNFSF18 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	<p>To detect hAITRL by sandwich ELISA (usin...</p>
Tested Applications	ELISA, WB	
Immunogen	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hAITRL. Human AITRL specific antibody was purified by affinity chromatography employing immobilized hAITRL matrix.	
Target	TNFSF18	
Form/Appearance	Lyophilized	<p>To detect hAITRL by Western Blot analysi...</p>
Concentration	batch dependent	
Storage	AITRL antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	<p>ELISA:Indirect:To detect hAITRL by indirect ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hAITRL.</p> <p>SandwichTo detect hAITRL by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with our biotinylated Anti-Human AITRL as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hAITRL.</p> <p>Western Blot:To detect hAITRL by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hAITRL is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.</p>	
Clonality	Polyclonal	
Uniprot ID	Q9UNG2	<p>To detect hAITRL by Western Blot analysi...</p>
NCBI	Q9UNG2	
Dilution Range	ELISA:Indirect:To detect hAITRL by indirect ELISA (using	