
Product Datasheet

DEFB103B Antibody (Biotin) (orb1272730)

Description

DEFB103B Antibody (Biotin)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Biotin

Tested Applications

ELISA, WB

Immunogen

Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hBD-3. Human BD-3 specific antibody was purified by affinity chromatography and then biotinylated.

Target

DEFB103B

Form/Appearance

Lyophilized

Concentration

batch dependent

Storage

BD-3 antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

ELISA: Sandwich: To detect hBD-3 by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.25 - 1.0 µg/mL of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our polyclonal Anti-Human BD-3 as a capture antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hBD-3. Western Blot: To detect hBD-3 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hBD-3 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.

Clonality

Polyclonal

Uniprot ID
P81534
NCBI
P81534
Dilution Range

ELISA: Sandwich: To detect hBD-3 by

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
 Research Triangle Park
 Durham, North Carolina
 27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

the detection of at least 0.2 - 0.4 ng/well of