

Product Datasheet

SARS-CoV-2 (COVID-19) S2 antibody [AbH201] (orb1272797)

Description	SARS-CoV-2 (COVID-19) S2 antibody [AbH201]
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 S2 protein
Preservatives	PBS,pH7.5
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Application notes	WB:1:1000-2000ELISA:1:5000-1:10000
Isotype	Human IgG
Clonality	Monoclonal
Dilution Range	WB:1:1000-2000ELISA:1:5000-1:10000
Expiration Date	12 months from date of receipt.