

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (orb1272816)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Preservatives	PBS,pH7.5
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Application notes	WB:1:1000-2000ELISA:1:5000-1:10000
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:1000-2000ELISA:1:5000-1:10000
Expiration Date	12 months from date of receipt.