
Product Datasheet

Ly6g Antibody (Biotin) (orb1272913)

Description

Ly6g Antibody (Biotin)

Species/Host

Rat

Reactivity

Mouse

Conjugation

Biotin

Tested Applications

FA, FC, ICC, IHC

Immunogen

Purified mouse BALB/c neutrophils.

Target

Ly6g

Preservatives

Liquid. In PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

0.5 mg/ml

Storage

Stable for at least 1 year after receipt when stored at -20°C.

Note

For research use only

Application notes

Flow Cytometry:(1:100).
 Immunohistochemistry:(1:200; frozen section / paraffin section) . Immunocytochemistry: (1:200). Functional (Blocking): 250ug / mouse (Optimal reagent to deplete neutrophils in vivo).

Isotype

IgG2b

Clonality

Monoclonal

Uniprot ID
[P35461](#)
NCBI
[547876](#)
Dilution Range

Flow Cytometry:(1:100).
 Immunohistochemistry:(1:200; frozen section / paraffin section) . Immunocytochemistry: (1:200). Functional (Blocking): 250ug / mouse (Optimal reagent to deplete neutrophils in vivo).

Expiration Date

12 months from date of receipt.