

Product Datasheet

NS5B Antibody [5B-12B7] (orb1272970)

Description	NS5B Antibody [5B-12B7]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FA, ICC, IP
Immunogen	Recombinant hepatitis C virus (HCV) NS5B RNA-dependent RNA polymerase.
Preservatives	Liquid. In PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Stable for at least 1 year after receipt when stored at -20°C.
Note	For research use only
Application notes	Western Blot: (1-2ug/ml) . Optimal conditions must be determined individually for each application.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	Western Blot: (1-2ug/ml) . Optimal conditions must be determined individually for each application.
Expiration Date	12 months from date of receipt.