

Product Datasheet

NS3 Antibody [1B6] (orb1272972)

Description	NS3 Antibody [1B6]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant hepatitis C virus (HCV) NS3 serine protease domain.
Preservatives	Liquid. In PBS containing 10% glycerol and 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Stable for at least 1 year after receipt when stored at -20°C.
Note	For research use only
Application notes	Immunocytochemistry: (paraformaldehyde fixation) (1:50). Western Blot: (1:500).
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Immunocytochemistry: (paraformaldehyde fixation) (1:50). Western Blot: (1:500).
Expiration Date	12 months from date of receipt.