

## Product Datasheet

### IL17A Antibody (FITC) (orb1272980)

<b>Description</b>	IL17A Antibody (FITC)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human IL-17A.
<b>Target</b>	IL17A
<b>Preservatives</b>	Liquid. 0.2um-filtered solution in PBS, pH 7.4. Contains no preservatives.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Stable for at least 1 year after receipt when stored at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Flow Cytometry: (1:20-1:500). Optimal conditions must be determined individually for each application.
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q16552</a>
<b>NCBI</b>	<a href="#">4504651</a>
<b>Dilution Range</b>	Flow Cytometry: (1:20-1:500). Optimal conditions must be determined individually for each application.
<b>Expiration Date</b>	12 months from date of receipt.