
Product Datasheet

**SARS-CoV-2 (COVID-19) Nucleocapsid Antibody
[3869] (orb1273018)**

Description

SARS-CoV-2 (COVID-19)

Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Coronavirus Nucleocapsid Antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Coronavirus Nucleocapsid Antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3869-3868 / 3869-3866
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3869-3868 / 3869-3866
Expiration Date	12 months from date of receipt.