

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3878] (orb1273024)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3878]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Coronavirus Nucleoprotein Antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Coronavirus Nucleoprotein Antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.