



## Product Datasheet

# SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3862] (orb1273031)

<b>Description</b>	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3862]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF
<b>Preservatives</b>	SARS-CoV-2 (COVID-19, 2019-nCoV) Nucleoprotein antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	SARS-CoV-2 COVID-19, 2019-nCoV) Nucleoprotein antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot Validation with Human Recom...



ELISA Test Antibodies: SARS-CoV-2 Nucleo...