



Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3851] (orb1273033)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3851]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Preservatives	SARS-CoV-2 (COVID-19, 2019-nCoV) Nucleoprotein antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	SARS-CoV-2 COVID-19, 2019-nCoV) Nucleoprotein antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Western Blot Validation with Human Recom...



ELISA Test. Antibodies: SARS-CoV-2 Nucle...