

Product Datasheet

Dengue Virus Antibody [5161] (orb1273034)

Description	Dengue Virus Antibody [5161]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Preservatives	Dengue Virus antibody is in a is formulated in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. Each vial contains 1.0 mlNo stabilizing proteins
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Dengue Virus MonoclonalDengue Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Suggested starting ranges are 1:10-1:50 for IFA, blotting and immunohistochemistry, and 1:20-1:200 for ELISA.
Clonality	Monoclonal
Dilution Range	Suggested starting ranges are 1:10-1:50 for IFA, blotting and immunohistochemistry, and 1:20-1:200 for ELISA.
Expiration Date	12 months from date of receipt.