

Product Datasheet

Parvovirus B19 Antibody [5191] (orb1273035)

Description	Parvovirus B19 Antibody [5191]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Parvovirus B19 antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Parvovirus B19 antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Suggested starting ranges are 1:10-1:50 for IF and 1:20-1:200 for ELISA
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Suggested starting ranges are 1:10-1:50 for IF and 1:20-1:200 for ELISA
Expiration Date	12 months from date of receipt.