

Product Datasheet

Mumps Virus Antibody [6007] (orb1273037)

Description	Mumps Virus Antibody [6007]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Mumps Virus antibody is in a is formulated in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. Each vial contains 1.0 ml. No stabilizing proteins
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Mumps Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Suggested starting ranges are 1:10-1:50 for Immunofluorescence and 1:20-1:200 for ELISA
Isotype	Mouse IgG2b
Clonality	Monoclonal
Dilution Range	Suggested starting ranges are 1:10-1:50 for Immunofluorescence and 1:20-1:200 for ELISA
Expiration Date	12 months from date of receipt.