

Product Datasheet

Rhinovirus Antibody [4967] (orb1273039)

Description	Rhinovirus Antibody [4967]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Rhinovirus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Rhinovirus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Rhinovirus Antibody [4965] can be used at 1:10-1:50 for Immunofluorescence. For ELISA use at 1:20-1:200.
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Rhinovirus Antibody [4965] can be used at 1:10-1:50 for Immunofluorescence. For ELISA use at 1:20-1:200.
Expiration Date	12 months from date of receipt.