

## Product Datasheet

### Canine Parvovirus Antibody [7351] (orb1273066)

|                            |                                                                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Canine Parvovirus Antibody [7351]                                                                                                                              |
| <b>Species/Host</b>        | Mouse                                                                                                                                                          |
| <b>Reactivity</b>          | Virus                                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                                                                                          |
| <b>Preservatives</b>       | Canine Parvovirus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added. |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                         |
| <b>Concentration</b>       | batch dependent                                                                                                                                                |
| <b>Storage</b>             | Canine Parvovirus antibody can be stored at -20°C, stable for one year.                                                                                        |
| <b>Note</b>                | For research use only                                                                                                                                          |
| <b>Isotype</b>             | Mouse IgG2a                                                                                                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                |