

Product Datasheet

Yellow Fever virus Antibody [3576] (orb1273114)

Description	Yellow Fever virus Antibody [3576]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Yellow Fever virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Yellow Fever virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.