

Product Datasheet

Herpes Virus (type 6) Antibody [2002] (orb1273122)

Description	Herpes Virus (type 6) Antibody [2002]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Preservatives	Herpes Virus antibody (type 6) is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Herpes Virus (type 6) antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.