

Product Datasheet

HIV Antibody [1981] (orb1273124)

Description	HIV Antibody [1981]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Preservatives	HIV antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	HIV antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.