

Product Datasheet

HIV Antibody [1943] (FITC) (orb1273129)

Description	HIV Antibody [1943] (FITC)
Species/Host	Mouse
Reactivity	Virus
Conjugation	FITC
Tested Applications	ELISA, IF, IHC, WB
Preservatives	HIV antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	HIV antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.