

Product Datasheet

Hepatitis B Virus Antibody [1892] (orb1273134)

Description	Hepatitis B Virus Antibody [1892]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Hepatitis B Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Hepatitis B Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.