

Product Datasheet

Hepatitis C Virus Antibody [1868] (orb1273141)

Description	Hepatitis C Virus Antibody [1868]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Preservatives	Hepatitis C Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Hepatitis C Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.