

Product Datasheet

Influenza A virus Antibody [1381] (orb1273173)

Description	Influenza A virus Antibody [1381]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Influenza A virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Influenza A virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.