

Product Datasheet

Parainfluenza Virus Antibody [1226] (orb1273188)

Description	Parainfluenza Virus Antibody [1226]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	Parainfluenza Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Parainfluenza Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.