

## Product Datasheet

# Respiratory Syncytial Virus Antibody [0681] (orb1273209)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Respiratory Syncytial Virus Antibody [0681]  |
| <b>Species/Host</b>        | Mouse  |
| <b>Reactivity</b>          | Virus  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, IF, NeA   |
| <b>Preservatives</b>       | Respiratory Syncytial Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added. |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | batch dependent  |
| <b>Storage</b>             | Respiratory Syncytial Virus antibody can be stored at -20°C, stable for one year.  |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | Mouse IgG2b  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |