

## Product Datasheet

### Epstein-Barr Virus Antibody [0256] (orb1273222)

<b>Description</b>	Epstein-Barr Virus Antibody [0256]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, IP, WB
<b>Preservatives</b>	Epstein-Barr Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Epstein-Barr Virus antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.