

Product Datasheet

Epstein-Barr Virus Antibody [0231] (orb1273224)

| | |
|----------------------------|---|
| Description | Epstein-Barr Virus Antibody [0231] |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | IF |
| Preservatives | Epstein-Barr Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added. |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Epstein-Barr Virus antibody can be stored at -20°C, stable for one year. |
| Note | For research use only |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |