

Product Datasheet

Herpes Simplex Virus Antibody [0113] (FITC) (orb1273228)

Description	Herpes Simplex Virus Antibody [0113] (FITC)
Species/Host	Mouse
Reactivity	Virus
Conjugation	FITC
Tested Applications	ELISA
Preservatives	Herpes Simplex Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Herpes Simplex Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.