

Product Datasheet

Vaccinia virus Antibody (HRP) (orb1273235)

Description	Vaccinia virus Antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Lister strain
Preservatives	Vaccinia virus antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Vaccinia virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.