

Product Datasheet

Respiratory Syncytial Virus Antibody (orb1273279)

Description	Respiratory Syncytial Virus Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human RSV isolate, confirmed
Preservatives	Respiratory Syncytial Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Respiratory Syncytial Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.