

Product Datasheet

Poliovirus Antibody (HRP) (orb1273285)

Description	Poliovirus Antibody (HRP)
Species/Host	Goat
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Pool of purified virus (types 1-3)
Preservatives	Poliovirus antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Poliovirus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.