

Product Datasheet

Poliovirus Antibody (FITC) (orb1273286)

Description	Poliovirus Antibody (FITC)
Species/Host	Goat
Reactivity	Virus
Conjugation	FITC
Tested Applications	WB
Immunogen	Pool of purified virus (types 1-3)
Preservatives	Poliovirus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Poliovirus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.