

Product Datasheet

Parainfluenza virus, type 3 Antibody (biotin) (orb1273289)

Description	Parainfluenza virus, type 3 Antibody (biotin)
Species/Host	Goat
Reactivity	Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Human isolate, type 3
Preservatives	Parainfluenza virus, type 3 antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Parainfluenza virus, type 3 antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.