

Product Datasheet

Parainfluenza virus, type 1 Antibody (orb1273295)

Description	Parainfluenza virus, type 1 Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Cantell strain
Preservatives	Parainfluenza virus, type 1 antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Parainfluenza virus, type 1 antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.