

Product Datasheet

Influenza A virus (H3N2) Antibody (FITC) (orb1273335)

Description	Influenza A virus (H3N2) Antibody (FITC)
Species/Host	Goat
Reactivity	Virus
Conjugation	FITC
Tested Applications	WB
Immunogen	Influenza A, strain Texas 1/77 (H3N2)
Preservatives	Influenza A virus (H3N2) antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Influenza A virus (H3N2) antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.