

Product Datasheet

Influenza A virus (H1N1) Antibody (biotin) (orb1273337)

Description	Influenza A virus (H1N1) Antibody (biotin)
Species/Host	Goat
Reactivity	Virus
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	Influenza A, strain USSR (H1N1)
Preservatives	Influenza A virus (H1N1) antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Influenza A virus (H1N1) antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.