

Product Datasheet

Influenza A virus (H1N1) Antibody (HRP) (orb1273338)

Description	Influenza A virus (H1N1) Antibody (HRP)
Species/Host	Goat
Reactivity	Virus
Conjugation	HRP
Tested Applications	IHC, WB
Immunogen	Influenza A, strain USSR (H1N1)
Preservatives	Influenza A virus (H1N1) antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Influenza A virus (H1N1) antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.