

Product Datasheet

HIV1 p24 Antibody (HRP) (orb1273346)

Description	HIV1 p24 Antibody (HRP)
Species/Host	Goat
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Purified native p24 from strain IIIB
Preservatives	HIV1 p24 antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	HIV1 p24 antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.