

## Product Datasheet

### HIV1 gp41 Antibody (HRP) (orb1273350)

<b>Description</b>	HIV1 gp41 Antibody (HRP)
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant (P. Pastoris) ectodomain of gp41 (glycosylated)
<b>Preservatives</b>	HIV1 gp41 antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	HIV1 gp41 antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.