

Product Datasheet

HIV1 gp120 Antibody (orb1273354)

Description	HIV1 gp120 Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified native gp120 from strain IIIB
Preservatives	HIV-1 gp120 antibody at 1 mg/ml in phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	HIV-1 gp120 antibody can be stored at 4°C or -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.