

Product Datasheet

HIV-1 (HTLV-3) Antibody (FITC) (orb1273357)

Description	HIV-1 (HTLV-3) Antibody (FITC)
Species/Host	Goat
Reactivity	Virus
Conjugation	FITC
Tested Applications	WB
Immunogen	Human isolate, confirmed
Preservatives	HIV-1 (HTLV-3) antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	HIV-1 (HTLV-3) antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.