

## Product Datasheet

### HIV-1 (HTLV-3) Antibody (orb1273358)

<b>Description</b>	HIV-1 (HTLV-3) Antibody
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Human isolate, confirmed
<b>Preservatives</b>	HIV-1 (HTLV-3) antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	HIV-1 (HTLV-3) antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.